



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-273


Anti NT5C

Product type	Primary antibodies
Immunogen	NT5C (AAH08183.1, 1 a.a. ~ 117 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MARSVRVLVDMGVLADFEAGLLRGFRRRFPEEPHVPLEQRRGFLAREQYRALRPDLADK VASVYEAPGFFLDLEIPGALDAVREMNDLPDLLKYHHCVGEKVWLPRPYSARGAA

Gene Alias DNT, DNT1, P5N2, PN-I, PN-II, UMPH2, cdN, dNT-1

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~13.30 kDa using NT5C transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of NT5C expression in transfected 293T cell line by NT5C rabbit polyclonal antibody

Lane 1: NT5C transfected lysate(13.30 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619