

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-355

Anti CYP46A1

Product type	Primary antibodies
Immunogen	CYP46A1 (NP_006659.1, 1 a.a. ~ 500 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence

MSPGLLLLLSAVLLAFGLCCTFVHRARSRYEHIPGPPRPSFLLGHLPCFWKKDEVGGRVL
QDVFLDWAKKYGPVVRVNVFHKTSVIVTSPESVKKFLMSTKYNKDSKMYRALQTVFGERL
FGQGLVSECNYERWHKQRRVIDLAFSRSSLVSLMETFNEKAEQLVEILEAKADGQTPVSM
QDMLTYTAMDILAKAAFGMETSMLLGAQKPLSQAVKLMLEGITASRNTLAKFLPGKRKQL
REVRESIRFLRQVGRDWVQRRREALKRGEVPADILTQILKAEEGAQDDEGLLDNFVTFF
IAGHETSANHLAFTVMELSRQPEIVARLQAEVDEVIGSKRYLDFEDLGRQLQYLSQVLKES
LRLYPPAWGTFRLLEEETLIDGVRVPGNTPLLSTYVMGRMDTYFEDPLTFNPDRFGPGA
PKPRFTYFPFSLGHRSCIGQQFAQMEVKVVMALLQRLEFRLVPGQRFGQLQEATLKPLD
PVLCTLRPRGWQPAPPPPPC


CP46, CYP46

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~56.80 kDa using CYP46A1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CYP46A1 expression in transfected 293T cell line by CYP46A1 rabbit polyclonal antibody.

Lane 1: CYP46A1 transfected lysate(56.80 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619