

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-380

Anti GMDS

Product type	Primary antibodies
Immunogen	GMDS (NP_001491.1, 1 a.a. ~ 372 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAHAPARCPSPARGSGDGEMGKPRNVALITGITGQDGSYLAEFLLKGYEVHGIVRRSSSF NTGRIEHLKYKNPQAHIEGNMKLHYGDLTDSTCLVKIINEVKPTEIYNLGAQSHVKISFDL AEYTADV DGVGTLRLLDVAVKTCGLINSVKFYQASTSELYGKVQEIPQKETTPFYPRSPYG AAKLYAYWIVVNFREAYNLFVNGILFNHESPRRGANFVTRKISRVAKIYLGQLECFSL GNLDAKRDWGHAKDYEAMWMLQNDEPEDFVIATGEVHSVREFVEKSFLHIGKTIVWEG KNENEVGRCKETGKVHVTVDLKYYRPTEVDLQGDCTKAKQKLNWKPRVAFDELVREMVH ADVELMRTNPNA


GMD, SDR3E1

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~41.90 kDa using GMDS transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of GMDS expression in transfected 293T cell line by GMDS rabbit polyclonal antibody.

Lane 1: GMDS transfected lysate(41.90 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of GMDS expression in mouse liver (Predicted Mw: 41.9500 kDa).

Western Blotting (Tissue lysate)



GMDS rabbit polyclonal antibody. Western Blot analysis of GMDS expression in mouse liver.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619