

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-390

Anti UBA2

Product type	Primary antibodies
Immunogen	UBA2 (NP_005490.1, 1 a.a. ~ 640 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence

MALSRGLPRELAEAVAGGRVLVVGAGGIGCELLKNLVLTGFSHIDLIDLTIDVSNLNKQ
 FLFQKKHVGRSKAQVAKESVLQFYFYPKANIVAYHDSIMNPDYNVEFFRQFILVMNALDNRA
 ARNHVNRMCLAADVPLIESGTAGYLGQVTTIKKGVTECYECHKPTQRTFPGCTIRNTPS
 EPIHCIVWAKYLFNQLFGEEDADQEVSPDRADPEAAWEPTAEARARASNEDGDIKRIST
 KEWAKSTGYDPVKLFTKLFKDDIRYLLTMDKLWRKRKPPVPLDWAQVQSQGEETNASDQQ
 NEPQLGLKDQQLDVKSYARLFSKSIETLRVHLAEKGDGAELIWDKDDPSAMDFVTSAN
 LRMHIFSMNMKSRFDIKSMAGNIIPAIATTNAVIAGLIVLEGLKILSGKIDQCRTIFLNK
 QPNPRKKLLVPCALDPPNPNCYVCASKPEVTVRLNVHKVTVLTLQDKIVKEKFAMVAPDV
 QIEDGKGITILISSEEGETEANNHKKLSEFGIRNGSRLQADDFLQDYTLINILHSEDLGK
 DVEFEVVGDAPEKVGPKQAEDAAKSITNGSDDGAQPSTSTAQEQDDVLIVDSDEEDSSNN
 ADVSEEERSRKRKLDEKENLSAKRSRIEQKEELDDVIALD


ARX, FLJ13058, HRIHFB2115, SAE2

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~71.20 kDa using UBA2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of UBA2 expression in transfected 293T cell line by UBA2 rabbit polyclonal antibody.

Lane 1: UBA2 transfected lysate(71.20 KDa).
 Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619