

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-396

Anti MRPL42

Product type	Primary antibodies
Immunogen	MRPL42 (NP_054769.1, 1 a.a. ~ 142 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAVAVKWVMSKRTILKHLFPVQNGALYCVCHKSTYSPLPDDYNCNVELALTSDGRTIVC YHPSVDIPYEHTKPIRPDPVHNNEETHDQVLKTRLEEKVEHLEEGPMIEQLSKMFFTTK HRWYPHGRYHRCRKNLNPPKDR


HSPC204, MRP-L31, MRPL31, MRPS32, PTD007, RPML31

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~16.70 kDa using MRPL42 transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of MRPL42 expression in transfected 293T cell line by MRPL42 rabbit polyclonal antibody.

Lane 1: MRPL42 transfected lysate(16.70 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619