

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-403

Anti DHRS

Product Primary antibodies

type

Immunogen DHRS11 (ENSP00000354182, 1 a.a. ~ 181 a.a) full-length human protein.

Form Liquid

Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MFSAIRSQHSGVDICINNAGLARPDTLTSGSTSGWKDMFNVNVLALSICTREAYQSMKER
NVDDGHIININSMSGHRVLPVTHFYATKYAVTALTEGLRQELREAQTHIRATCISPG
VVETQFAFKLHDKDPEKAAATYEQMKCLKPEDVAEAVIYVLSTPAHIQIGDIQMRPTEQV
T


ARPG836, FLJ39232, MGC4172, SDR24C1

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~20.00 kDa using DHRS11 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of DHRS11 expression in transfected 293T cell line by DHRS11 rabbit polyclonal antibody.

Lane 1: DHRS11 transfected lysate(20.00 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of DHRS11 expression in mouse intestine (Predicted Mw: 28.3080 kDa).

Western Blotting (Tissue lysate)



DHRS11 rabbit polyclonal antibody. Western Blot analysis of DHRS11 expression in mouse intestine.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619