

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-406

Anti AKR1B1

Product Primary antibodies

type

Immunogen AKR1B1 (AAH00260.1, 1 a.a. ~ 316 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MASRLLLNNGAKMPILGLGTWKSPPGQVTEAVKVAIDVGYRHIDCAHVYQNEVEVGVAIQ
EKLREQVVKREELFIVSKLWCTYHEKGLVKGACQKTLSDLKLDYLDLYLIHWPTGFKPGK
EFFPLDESGNVVPSDTNILDWAAMEELVDEGLVKAIGISNFNHLQVEMILNKPGKLYKP
AVNQIECHPYLTQEKLQYQCSKGIVVTAYSPLGSPDRPWAKPEDPSLLEDPRIKAIKAAK
HNKTTAQVLIRFPMQRNLVVIPKSVTPERIAENFKVDFELSSQDMTLLSYNRNWRVCA
LLSCTSHKDYPFHEEF


ADR, ALDR1, ALR2, AR, MGC1804

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~34.76 kDa using AKR1B1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of AKR1B1 expression in transfected 293T cell line by AKR1B1 rabbit polyclonal antibody.

Lane 1: AKR1B1 transfected lysate(34.76 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619