



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-473


Anti CDKN1A

Product type	Primary antibodies
Immunogen	CDKN1A (NP_000380.1, 1 a.a. ~ 164 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MSEPAGDVRQNPCGSKACRRLFGPVDSEQLSRDCDALMAGCIQEARERWNFDFVTETPLE GDFAWERVRLGLPKLYLPTGPRRGRDELGGRRPGTSPALLQGTAEDHVDLSLSCTLV PRSGEQAEGSPGGPGDSQGRKRRQTSMTDFYHSKRRLIFSKRKP

Gene Alias CAP20, CDKN1, CIP1, MDA-6, P21, SDI1, WAF1, p21CIP1

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~18.04 kDa using CDKN1A transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CDKN1A expression in transfected 293T cell line by CDKN1A rabbit polyclonal antibody.

Lane 1: CDKN1A transfected lysate(18.04 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619