

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-521

Anti YWHAG

Product type Primary antibodies

Immunogen YWHAG (NP_036611.2, 1 a.a. ~ 247 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MVDREQLVQKARLAEQAERYDDMAAAMKNVTELNEPLSNEERNLLSVAYKNVVGARRSSW
RVISSIEQKTSADGNEKKIEMVRAYREKIEKELEAVCQDVLSDNYLIKNCSETQYESK
VFYLLMKGDYYRYLAEVATGEKRATVVESSEKAYSEAHEISKEHMQPTHPIRLGLALNYS
VFYYEIQNAPEQACHLAKTAFDDAIAELDTLNEDSYKDSTLIMQLLRDNLTLWTSDQQDD
DGEGNN


Gene Alias 14-3-3GAMMA

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~28.30 kDa using YWHAG transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of YWHAG expression in transfected 293T cell line by YWHAG rabbit polyclonal antibody.

Lane 1: YWHAG transfected lysate(28.30 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619