



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-536

Anti CDC42EP2

Product type	Primary antibodies
Immunogen	CDC42EP2 (NP_006770.1, 1 a.a. ~ 210 a.a) full-length human protein
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MSTKVPIYLKRGSRKGKKEKLRDLLSSDMISPLLGDFRHTIHIGSGGGSDMFGDISFLQG KFHLLPGTMVEGPEEDGTFDLPFQFTRTATVCGRELPDGSPLLKNAISLPVIGGPQALT LPTAQAPPKPPRLHLETPQSPQEGGSVDIWRI PETGSPNSGLTPESGAEEPFLSNASSL LSLHVDLGPSILDDVLQIMDQDLDSMQIPT


Gene Alias BORG1, CEP2

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~22.50 kDa using CDC42EP2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CDC42EP2 expression in transfected 293T cell line by CDC42EP2 rabbit polyclonal antibody.

Lane 1: CDC42EP2 transfected lysate(22.50 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619