



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-537

Anti CDC42EP3

Product type	Primary antibodies
Immunogen	CDC42EP3 (NP_006440.2, 1 a.a. ~ 254 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MPAKTPIYLKAANNKKGKKFKLRDILSPDMISPPLGDFRHTIHIGKEGQHDVFGDISFLQ GNYELLPGNQEK AHLGQFPGHNEFFRANSTSDSVFTETPSPVLKNAISLPTIGGSQALML PLLSPVTFNSKQESFGPAKLPRLSCEPVMEEKAQEKSSLLENGTVHQGDTSWGSSGSASQ SSQGRDSHSSSLSEQYPDWPAEDMFDHPTPCELIKGKTKSEESLSDLTGSLLSLQLDLGP SLLDEV LNVMDKNK


Gene Alias BORG2, CEP3, FLJ46903, UB1

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~27.70 kDa using CDC42EP3 transfected 293T lysate

Western Blotting (Transfected lysate)




Western Blot analysis of CDC42EP3 expression in transfected 293T cell line by CDC42EP3 rabbit polyclonal antibody.

Lane 1: CDC42EP3 transfected lysate(27.70 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of CDC42EP3 expression in A-431 (Predicted Mw: 27.6780 kDa).

Western Blotting (Cell lysate)



CDC42EP3 rabbit polyclonal antibody. Western Blot analysis of CDC42EP3 expression in A-431.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619