

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-590

Anti POU2F1

Product type Primary antibodies

Immunogen POU2F1 (AAH03571.1, 1 a.a. ~ 743 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MNNPSETSKPSMESGDGNTGTQTNGLDFQKQPVPVGGAISTAQAQAFGLGHLHQVQLAGTS
LQAAAQSLNVQSKSNEESGDSQQPSQPSQQPSVQAAPQTQLMLAGGQITGLTLTPAQQQ
LLLQQAQAQAQLLAAAVQQHSASQQHSAAGATISASAATPMTQIPLSQPIQIAQDLQQLQ
QLQQQNLNLQQFVLVHPTTNLQPAQFIISQTPQGQQGLLQAQNLTLQLPQQSQANLLQSQ
PSITLTSQPATPTRTIAATPIQTLPQSQSTPKRIDTPSLEEPSDLEELEQFAKTFKQRRRI
KLGFTQGQDVGLAMGKLYGNDFSQTTISRFEALNLSFKNMCKLKPILLEKWLND AENLSSDS
SLSSPSALNSPGIEGLSRRRKRTSIETNIRVALEKSFLNQKPTSEEITMIADQLNMEK
EVIRVWFCNRRQKEKRINPPSSGGTSSSPIKAIFPSPTSLVATTPSLVTSSAATTLTVSP
VLPLTSAAVTNLSVTGTSDTTSNNTATVISTAPPASSAVTSPSLSPSPSASASTSEASSA
SETSTTQTTSTPLSSPLGTSQVMVTASGLQTA AAAALQGAAQLPANASLAAMAAAAGLNP
SLMAPSQFAAGGALLSLNPGTSLGALSPALMSNSTLATIQALASGGSLPITSLDATGNLV
FANAGGAPNIVTAPLFLNPQNLSLLTSNPVSLVSA A AASAGNSAPVASLHATSTSAESIQ
NSLFTVASASGAASTTTTASQAQ


Gene Alias OCT1, OTF1

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~76.50 kDa using POU2F1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of POU2F1 expression in transfected 293T cell line by POU2F1 rabbit polyclonal antibody.

Lane 1: POU2F1 transfected lysate(76.50 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C or one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619