

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-WDR20-117

Anti WDR20

Product type	Primary antibodies
Immunogen	WDR20 (AAH28387.1, 1 a.a. ~ 569 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MATEGGGKEMNEIKTQFTTREGLYKLLPHSEYSRPNRVPFNSQGSNPVRVSFVNLNDQSG NGDRLCFNVGRELYFYIYKGVKAAADLSKPIDKRIYKGTQPTCHDFNHLTATAESVLLV GFSAGQVQLIDPIKKETSKLFNEERLIDKSRVTCVKWVHGSESLFLVAHSSGNMYLYNVE HTCGTTAPHYQLLKQGESFAVHTCKSKSTRNPLLKWTVGEGALNEFAFSPDGKFLACVSQ DGFLRVFNFDSELHGTMKSYFGLLCVCWSPDGKYIVTGGEDDLTVVWSFVDCRVIARG HGHKSWVSVVAFDPYTTVEEGDPMFSGSDEDFQDLLHFGRDRANSTQSRLSKRNSTD RPVSVTYRFGSVGQDTQLCLWDLTEDILFPHQPLSRARTHNTVMNATSPAGSNGNSVTT PGNSVPPPLPRNSLPHSAVSNAGSKSSVMDGAIASGVSKFATLSLHDKERHHEKDHKR NHSMGHISSKSSDKLNLVTKTCTDPAKTLGTPLCPRMEDVPLLEPLICKKIAHERLTVLI FLEDCIVTACQEGFICTWGRPGKVVSFNP
Gene Alias	DMR, FLJ33659, MGC33177, MGC33183
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~62.90 kDa using WDR20 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of WDR20 expression in transfected 293T cell line by WDR20 rabbit polyclonal antibody.

Lane 1: WDR20 transfected lysate(62.90 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619