

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.


Catalog No. CNP-YARS2-284

Anti YARS2

Product type	Primary antibodies
Immunogen	YARS2 (NP_001035526.1, 1 a.a. ~ 477 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAAPILRSFSWGRWSGTLNLSVLLPLGLRKAHSGAQGLLAAQKARGLFKDFFPETGTKIE LPELFDRTGASFPQTIYCGFDPTADSLHVGHLLALLGLFHLQRAGHNIALVGGATARLG DPSGRTERKEREALETERVRANARALRLGLEALAAHQQLFTDGRSWGSSFTVLDNSAWYQKQ HLVDFLAAVGGHFRMGTLRSRQSVQLRLKSPEGMSLAEFFYQVLQAYDFYYLFQRYGCRV QLGGSDQLGNIMSGYEFINKLTGEDVFGITVPLITSTTGAKLGKSAGNAVWLNDRDKTSPF ELYQFFVRQPDDSVRYLKLFTFLPLPEIDHIMQLHVKEPERRGPQKRLAAEVTCLVHGR EGLDSAKRCTQALYHSSIDALEVMDSQELKELFKEAPFSEFFLDPGTSVLDTCRKANAIP DGPRGYRMITEGGVSINHQVTNPESVLIVGQHILKNGLSLLKIGKRNFYIWKWLQL

Gene Alias	CGI-04, FLJ13995, TYRRS, mt-TyrRS
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~53.20 kDa using YARS2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of YARS2 expression in transfected 293T cell line by YARS2 rabbit polyclonal antibody.

Lane 1: YARS2 transfected lysate(53.20 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619