

MONOCLONAL ANTIBODY

For research use only, Not for diagnostic use.

Catalog No. KYU-HT-M005

Anti Plumbagin

[PL; 5-hydroxy-2-methyl-1,4-naphthoquinone]

BACKGROUND

Plumbagin is a naphthoquinone isolated from Plumbago zeylanica.

Product type	Primary antibody
Immunogen	PL(Plumbagin)-BSA
Raised in	Mouse (Balb/c)
Myeloma	SP2/0
Clone number	3A3
Isotype	IgG1 k
Source	Serum free culture supernatant
Purification	Protein G
Buffer	PBS containing 0.02% NaN ₃ as a preservative
Concentration	1mg / mL
Volume	100 UL
Label	Unlabeled
Specificity	plumbagin
Cross reactivity	-
Storage	Store below -20 degrees (below -70 degrees for prolonged storage) Aliquot to avoid cycles of freeze / thaw.

Application notes • **ELISA** (Ref.1)

Recommended dilutions
Other applications have not been tested.
Optimal dilutions/concentrations should be determined by the end user.

References
1) Sakamoto et al., Development of an enzyme-linked immunosorbent assay (ELISA) using highly-specific monoclonal antibodies against plumbagin, Anal Chim Acta. 2008 Jan 21;607(1):100-5. PMID: [18155415](#)

For research use only, Not for diagnostic use.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619