

**MONOCLONAL ANTIBODY**

*For research use only. Not for clinical diagnosis.*

**Catalog No. MKM-M01**

## Anti basic FGF, clone bFM-1

<b>Product type</b>	Primary antibodies
<b>Host</b>	Mouse
<b>Source</b>	Ascites (Affinity purified)
<b>Form</b>	Liquid
	PBS (pH 7.4)
<b>Volume</b>	100 µl
<b>Concentration</b>	1.0 mg/ml
<b>Specificity</b>	Basic FGF
<b>Antigen</b>	Purified bovine basic FGF
<b>Clone</b>	bFM-1
<b>Isotype</b>	IgG1

**Application notes** Radioimmunoassay, Cell biology: neutralizing, not recommended for western blotting  
**Recommended use**

**Recommended dilutions**

Optimal dilutions/concentrations should be determined by the end user.

**Staining Pattern**

**Cross reactivity** Human, Bovine, Mouse, Rat

**Storage** Store below -20°C (below -70°C for prolonged storage).  
Aliquot to avoid cycles of freeze/thaw.

**References** 1) Matsuzaki K et al. Monoclonal antibodies against heparin-binding growth factor II/basic growth factor that block its biological activity: invalidity of the antibodies for tumor angiogenesis. P.N.A.S. 86: 9911-9915 (1989)



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619