

MONOCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. MKM-M09

Anti SnoN, clone 3'-1

Product type	Primary antibodies
Host	Mouse
Source	Ascites (Affinity purified)
Form	Liquid
	PBS (pH 7.4)
Volume	100 µl
Concentration	1.0 mg/ml
Specificity	SnoN
Antigen	Recombinant human SnoN
Clone	3'-1
Isotype	IgG1 with λ-type light chains

Application notes WB, ELISA, IP
Recommended use

Recommended dilutions

Western blotting: 1/100

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity Human

Storage Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

References

- 1) Shinagawa, T. et al. The sno gene, which encodes a component of the histone deacetylase complex, acts as a tumor suppressor in mice. EMBO J. 19: 2280-2291 (2000)
- 2) Kokura K. et al. The Ski-binding protein C184M negatively regulate TGF-β signaling by sequestering the Smad proteins to cytoplasm. J. Biol. Chem. 278: 20133-20139 (2003)



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619