

MONOCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. MKM-M14

Anti PF-4, clone 2D6B

Product type	Primary antibodies
Host	Mouse
Source	Ascites (Affinity purified)
Form	Liquid
	PBS (pH 7.4)
Volume	100 µl
Concentration	1.0 mg/ml
Specificity	PF-4
Antigen	Human platelet factor 4 (PF-4)
Clone	2D6B
Isotype	IgG1 with κ-type light chains

Application notes WB, IP

Recommended use

Recommended dilutions

Western blotting: 1/100

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity Not examined

Storage Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

References



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619