



MONOCLONAL ANTIBODY

Catalog No. PUH-031N-A

Anti human TNF- α

Product type	Primary antibodies
Host	Mouse
Source	Culture supernatant
Form	
Volume	100ul
Concentration	1 mg/ml
Specificity	human TNF- α
Antigen	Human recombinant TNF- α
Clone	ETNF-10
Isotype	IgG1k

Application notes

WB

Recommended use

Recommended dilutions

Western Blot: 1/100

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Human

Storage

Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

References

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti human TNF- α



←17kDa



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619