

MONOCLONAL ANTIBODY

*For research use only, Not for diagnostic use.*

**Catalog No. CAC-SHI-EXO-M02-EX**

## Anti CD63 [ Clone : 8A12]

### BACKGROUND

CD63 (also known as LAMP-3, Melanoma-associated antigen ME491, TSPAN30, MLA1 and OMA81H) is a cell surface glycoprotein which belongs to the tetraspanin superfamily. CD63 is known to complex with integrins. CD63 is expressed on activated platelets, monocytes and macrophages, and is weakly expressed on granulocytes, T cell and B cells.

<b>Product type</b>	Primary antibody
<b>Immunogen</b>	Human CD63(a.a 104-202)
<b>Rased in</b>	Mouse
<b>Myeloma</b>	-
<b>Clone number</b>	8A12
<b>Isotype</b>	IgG2a
<b>Source</b>	Culture supernatant
<b>Purification</b>	Affinity purified by Protein A
<b>Buffer</b>	Phosphate buffered saline containing 0.1% NaN <sub>3</sub> as a preservative
<b>Concentration</b>	1 mg / mL
<b>Volume</b>	100 ul
<b>Label</b>	Unlabeled
<b>Specificity</b>	Human CD63
<b>Cross reactivity</b>	Human. Other species is not tested.
<b>Storage</b>	Store below - 40°C Aliquot to avoid cycles of freeze/thaw..


---

<b>Application notes</b>	• <b>Western blotting:</b> 1/500 - 1/2,000 (0.5 - 2ug/mL, non-reductive condition)
<b>Recommended dilutions</b>	• <b>Immunoprecipitation:</b> 1/200 – 1/1,000 (1- 5ug/mL) * Experiment example : Serum and Culture supernatant Exosome


Other applications have not been tested.  
Optimal dilutions/concentrations should be determined by the end user.

<b>References</b>	-
-------------------	---

## ANTIBODY CHARACTERIZATION




**Fig.1 Western blot** analysis of CD63 in C32 cell Exosome with anti-CD63 (8A12) antibody




**Fig.2 Immunoprecipitation** in Human serum exosome

### Protocol

Sample	
Serum	150uL
+	
Anti-CD63(8A12) 1ug (Solid phase Magnetic bead)	
PBS	150uL





**Fig.3 Immunoprecipitation** in C32 cell (Melanoma cell line) Exosome

**Protocol**

Sample	
C32 cell Exosome	2 ug
+	
Anti-CD63(8A12)	1ug (Solid phase Magnetic bead)
PBS+1%BSA	200uL(Total)



Rotated at 4°C o/n



Washed

SDS-PAGE sample buffer 15uL



Denatured by heating

SDS-PAGE (3uL)  
WB (Anti-CD63(8A12) : 1ug/mL)

**RELATED PRODUCT:**

Product Name	Application	Quantity	Cat#
Anti CD9 (12A12) Monoclonal Antibody	WB / IP / ELISA	100 ug	CAC-SHI-EXO-M01-EX
Anti CD63 (8A12) Monoclonal Antibody	WB / IP	100 ug	CAC-SHI-EXO-M02-EX
Anti CD81 (12C4) Monoclonal Antibody	WB / IP	100 ug	CAC-SHI-EXO-M03-EX

For research use only, Not for diagnostic use.

Manufactured by  **SHIONOGI & CO., LTD.**



**COSMO BIO Co., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619