



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. SK-T01-005

Anti-human metastin antiserum

Product type	Primary antibodies
Host	Rabbit
Source	Serum
Form	Liquid
	0.1% NaN ₃ as a preservative.
Volume	100 µl
Concentration	
Specificity	
Antigen	Full-length Human metastin -GST fusion protein
Isotype	IgG

Application notes WB, IHC, IP, ELISA
Recommended use

Recommended dilutions

Western blotting: 1/500 to 1/1,000

Predicted molecular weight: 60 kDa

Immunohistochemistry: 1 µl

Immunoprecipitation: 1/100 to 1/500

ELISA: 1/2,000

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity Not tested


Storage Store below -20°C (below -70°C for prolonged storage).
Aliquot to avoid cycles of freeze/thaw.

References



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti-human metastin antiserum



For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619