

POLYCLONAL ANTIBODY

Catalog No. YCU-PS-M3

Anti muscarinic acetylcholine receptor m3

Product type	Primary antibodies
Host	Rabbit
Source	Whole Serum
Form	Liquid
Volume	1 ml
Concentration	
Specificity	Muscarinic acetylcholine receptor m3
Antigen	Muscarinic acetylcholine receptor m3 (S302-Q473)
Isotype	

Application notes IP Not yet tested in other applications.

Recommended use

Recommended dilutions

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity Human, Rat

Storage Store below -20°C (below -70°C for prolonged storage). Aliquot to avoid cycles of freeze/thaw.

References

- 1) K. Shiozaki, E. Iseki, H. Uchiyama, Y. Watanabe, T. Haga, K. Kameyama, T. Ikeda, T. Yamamoto, K. Kosaka Alterations of muscarinic acetylcholine receptor subtypes in diffuse Lewy body disease: relation to Alzheimer's disease J Neurol Neurosurg Psychiatry. 1999 Aug;67(2):209-13.
- 2) Kazumasa Shiozaki, Eizo Iseki, Hiroaki Hino, Kenji Kosaka Distribution of m1 muscarinic acetylcholine receptors in the hippocampus of patients with Alzheimer's disease and dementia with Lewy bodies-an immunohistochemical study J Neurol Sci. 2001 Dec 15;193(1):23-8.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619