



COSMO BIO CO., LTD.

Inspiration for Life Science

Anti ovomucoid monoclonal antibody (CB-2)

Product No.:	CBN-CH-006-EX
Origin:	Mouse
Antigen:	ovomucoid
Subclass:	IgG1
Quantity:	0.1 mg (2.0 mg / 1mL)
Product form:	protein A purified IgG in PBS,0.095%NaN ₃
Storage:	Store below -20°C
Culture Medium:	DMEM,10% Calf Serum,
Myeloma cell:	P3X63Ag8U1
Application:	ELISA Working dilution is 2,000 times The maximum positivity dilution is 300,000 times

FOR RESEARCH LABORATORY USE ONLY



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN

TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp