



COSMO BIO CO., LTD.

Inspiration for Life Science

Anti gliadin monoclonal antibody(CE-3)

Product No.:	CBN-CH-013-EX
Origin:	Mouse
Antigen:	Wheat gliadin
Subclass:	IgG1
Quantity:	0.1 mg (2.0 mg / 1mL)
Product form:	protein A purified IgG in PBS,0.095%NaN ₃
Storage:	Store below -20°C
Culture Medium:	RPMI1640,10% Calf Searum, 55uM mercaptoethanol 1mM Sodium Pyruvate, 1×non-essential amino acid
Myeloma cell:	P3X63Ag8U1
Application:	Western Blot Working dilution is 7,800 times ELISA Working dilution is 7,800 times The maximum positivity dilution is 1,000,000 times

FOR RESEARCH LABORATORY USE ONLY



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN

TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp