



COSMO BIO Co., LTD.

Inspiration for Life Science

Recombinant protein

For research use only, Not for diagnostic use

Catalog No. KN-TOYU-M01

Transthyretin(His-Tag)

Product type	Recombinant Protein / Transthyretin (His-Tag)
Sequence	MRGSHHHHHHSGSGPTGTGESKCPLMKVLDVAVRGSPAINVAVHVFRKAADDTWEPFASG KTSESGELHGLTTEEEFVEGIYKVEIDTKSYWKALGISPFHEHAEEVFTANDSGPRRYTIA ALLSPYSYSTTAVVTNPKE
Source	E.coli
Form	Lyophilized
Volume	1 mg / tube
Storage Conditions	Store at 4°C. After reconstituted, store no longer than 1-2 days at 4°C, unless preservative is added. For long term storage, aliquot and store reconstituted product frozen at -20°C or lower. Aliquot to avoid cycles of freeze/thaw.

For research use only, Not for diagnostic use.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9613

FAX : +81-3-5632-9619