

Recombinant bTIMP-1

This is a purified recombinant bovine TIMP-1 (r-bTIMP-1) from CHO cells transfected with bTIMP-1 cDNA.

Concentration	100 ug/mL
Volume	1 mL
Buffer	PBS(-): 10mM Sodium phosphate buffer, pH7.3 containing 0.15M NaCl
Storage	-20 centigrade
Product No.	Y-10
Reference	1)Freudenstein, J. et al. mRNA of bovine tissue inhibitor of metalloproteinase: Sequence and expression in bovine obarian tissue. Biochem. Biophys. Res. Commun. 171, 250-256, 1990.

Note For laboratory use only.

Revised in Jan., 1998.
Revised in Oct., 2001.

Distributor

COSMO BIO

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
TEL.: +81-3-5632-9605
FAX.: +81-3-5632-9614
URL: <http://www.cosmobio.co.jp>
e-mail: export@cosmobio.co.jp

Manufacture



530 Chokeiji Takaoka Toyama 933-8511, Japan
TEL:81-766-26-4423 / FAX:81-766-26-4437
e-mail: msts@daiichi-fcj.co.jp
Web URL: <http://www.daiichi-fcj.co.jp>