

## Recombinant hTIMP-2

This is a purified recombinant human TIMP-2 (r-hTIMP-2) from CHO cells transfected with hTIMP-2 cDNA.

Concentration	100 ug/mL
Volume	1 mL
Buffer	PBS(-): 10mM Sodium phosphate buffer, pH7.3 containing 0.15M NaCl
Storage	-20 centigrade
Product No.	Y-9
Reference	1)Boone, T. C. et al. cDNA cloning and expression of a metalloproteinase inhibitor related to tissue inhibitor of metalloproteinases. Proc. Natl. Acad. Sci. USA, 87, 2800-2804, 1990.  2)Aoki, T. et al. Expression and large scale production of human tissue inhibitor of metalloproteinases-2 from recombinant chinese hamster ovary cells by high density continuous culture. Connect. Tissue, 26, 281-290, 1995. (In Japanese)  3)Kumagai, K. et al. Inhibition of matrix metalloproteinases prevents allergen-induced airway inflammation in a murine model of asthma. J. Immunol., 162, 4212-4219, 1999.

Note For laboratory use only.

Revised in June, 1999.  
Revised in Oct., 2001.

### Distributor

**COSMO BIO**

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,  
KOTO-KU. TOKYO 135-0016, JAPAN  
TEL.: +81-3-5632-9605  
FAX.: +81-3-5632-9614  
URL: <http://www.cosmobio.co.jp>  
e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

### Manufacture



530 Chokeiji Takaoka Toyama 933-8511, Japan  
TEL:81-766-26-4423 / FAX:81-766-26-4437  
e-mail: [msts@daiichi-fcj.co.jp](mailto:msts@daiichi-fcj.co.jp)  
Web URL: <http://www.daiichi-fcj.co.jp>