

Anti-Medaka Vitellogenin (F10-2)

Mouse monoclonal antibody

Background:	Vitellogenin, the egg yolk precursor protein, is induced by estrogens. On the other hand, environmental estrogens are known to disrupt endocrine system in animals. Recently, vitellogenin has proven to be an ideal marker for environmental estrogens.
Code:	AVMF10-2
Lot:	
lot size:	1 mL
Isotype:	IgG1 κ chain
Purity :	>80% by SDS-PAGE
Concentration:	0.1mg/ml
Source:	purified Medaka vitellogenin
Appearance:	Clear liquid
Storage:	Short term, 4° C Long term, -20° C Avoid freeze thaw cycles
Preservative:	0.1% Sodium azide
Presentation:	Supplied as a liquid in 0.01M Tris,0.15M NaCl,0.1% Sodium azide,pH:7.5

For research use only; not for use as a diagnostic.

COSMO BIO CO.,LTD

TOYO EKIMAE BLDG, 2-20,TOYO 2CHOME,KOTO-KU,TOKYO 135-0016,JAPAN
TEL: +81-3-5632-9605 FAX: +81-3-5632-9614 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp