

Anti-Fathead Minnow Vitellogenin : 883 (Mouse monoclonal antibody)

Code :	EBT-FHM883-EX
Lot :	
Lot size :	1 mL
Isotype :	IgG
Purity :	>80% by SDS-PAGE
Concentration :	0.1 mg /mL
Background:	Vitellogenin, the egg yolk precursor protein, is induced by estrogens. On the other hand, environmental estrogens are known to disrupt endocrine system in animals. Recently, vitellogenin has proven to be an ideal marker for environmental estrogens.
Source :	Purified plasma of 17 β -estradiol treated female fathead minnow
Appearance :	Clear liquid
Storage :	Short term, 4°C Long term, -20°C Avoid freeze thaw cycles
Preservative :	0.1% Sodium azide
Presentation :	Supplied as a liquid in 0.01M Tris, 0.15M NaCl, 0.1% Sodium azide, pH:7.5

Distributor



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

<http://www.cosmobio.co.jp>

Phone : +81-3-5632-9617

e-mail : export@cosmobio.co.jp

FAX : +81-3-5632-9618