



Progesterone-3(E)-CMO-AB-HRP

Product type

Form	Liquid
Volume	1 ml(1000Test)
Concentration	

Application notes

ELISA

Recommended use

Recommended dilutions

ELISA : 50 µl/Well use

Storage

Optimal dilutions/concentrations should be determined by the end user.

Maintain at 4°C-8°C (-20°C for long term strage) in undiluted aliquots. Avoid freeze/thaw cycles

References

* The Cross reactivities of Anti Progesterone-3(E)-CMO (FKA-302E) by ELISA with Progesterone-3(E)-CMO-AB-HRP (FKA-301)

Steroid

Cross Reaction

Progesterone	100 %
5Alpha-Pregnanedione	12.50 %
11Alpha-OH-Progesterone	5.30 %
20Alpha-OH-Progestronee	0.20 %
Pregnenolone	2.00 %
DOC	0.01 %
17Alpha-OH-Progesterone	0.01 %
Corticosterone	0.01 %
Cortisol	0.01 %
Aldosterone	0.01 %



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619