



Anti Cortisol-3 polyclonal antibody

Product No. FKA-404

The Cross reactivities of anti Cortisol-3 (FKA-404) by RI Assay

Antigen	Cortisol-3-CMO-BSA	
Antibody[FKA-404]	lyophilized Serum	25μl / Vial
Dilution of Serum	2000test for RI Assay	250μl / tube use
Binding of Total	35%	
Cortisol dose of B0/2	400 pg	
Strage/Handling	Maintain at 4°C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.	

Steroid	Cross Reaction
Cortisol	100 %
11-Deoxycortisole	16.6 %
Cortisone (Compd. E)	11.0 %
Corticosterone (Compd.E)	1.7 %
17alpha-Hydroxy-11-deoxy-progesterone (Compd. A)	0.9 %
17alpha-Hydroxy-orogesterone	0.9 %
The other steroids	0.0 %

For research use only; not for use as a diagnostic.

