



COSMO BIO Co., LTD.

Inspiration for Life Science

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. FKA-428

Anti Aldosterone-3-CMO-BSA

Product type	Primary antibodies
Host	Rabbit
Source	Serum
Form	Lyophilized powder
Volume	1 µl (2000 TEST)
Concentration	
Specificity	Aldosterone-3-CMO-BSA *
Antigen	Aldosterone-3-CMO-BSA
Isotype	

Application notes Radio Immunoassay

Recommended use

Recommended dilutions

RIA : 250 µl/tube use

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity Human, Equine, Bovine, Mammalian, Fish

Storage Maintain at 4°C-8°C (-20°C for long term storage) in undiluted aliquots. Avoid freeze/thaw cycles



References

*The Cross reactivities of anti Aldosterone-3-CMO-BSA (FKA-428) by RI Assay

Steroid	Cross Reaction
Aldosterone	100%
DOC	0.04%
Cortisone	0.01%
Cortisol	0.004%
Corticosterone	0.004%
Progesterone	0.07%
Testosterone	0.02%
Androstenedione	0.05%
Estradiol	0.002%

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: http://www.cosmobio.co.jp/index_e.asp e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619