



## PRODUCT DATA SHEET

**Product:** Anti-Amelogenin, Polyclonal

**Cat.No:** AB 01 ( 100  $\mu$  l/vial)

**Species Reactivity:**

Bovine ,rat, and mouse. Other species no tested .

**Immunogen**

Bovine amelogenin .

**Storage and Stability:**

Store at  $-20^{\circ}\text{C}$  or below.

Aliquot to avoid cycles of freeze/thaw

**Application and Suggested Dilutions:**

Immunohistochemistry:

staining with PAP or ABC methods . Suggested dilution: 100–500 fold .

Western Blot :

Use at a dilution of 1:1,000 to 1:10,000 .

The optimal dilution for a specific application should be determined by the researcher.

**Limitations:**

For in vitro research use only. Not for use in diagnostics or in humans.

**Distributor**



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME  
KOTO-KU, TOKYO 135-0016, JAPAN  
TEL : +81-3-5632-9617  
FAX : +81-3-5632-9618  
URL : <http://www.cosmobio.co.jp/>  
e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

**Manufacturer**



HOKUDO CO.,LTD

4-28,West10,Hachiken9jyo,Nishi-ku,  
Sapporo,063-0849,JAPAN  
TEL. : +81-11-641-7507  
FAX. : +81-11-644-9209  
URL : <http://www.hokudo.co.jp>  
E-mail : [info@hokudo.co.jp](mailto:info@hokudo.co.jp)