



PRODUCT DATA SHEET

Catalog Number: HKD-AB0217C5-EX

Product: Anti-Osteopontin

Clone: No. 17C5

Quantity: 50µl / vial

Host:

Mouse monoclonal antibody

Species Reactivity:

human, others not tested.

Immunogen:

human native osteopontin.

Format:

50ul of purified mouse monoclonal antibody by ammonium sulfate.

Buffer: PBS(-)

Storage and Stability:

Stable for 2 years when stored at less than -20°C.

Aliquot to avoid freeze/thaw cycles.

Application and Suggested Dilutions:

The optimal dilution for a specific application should be determined by the researcher.

For in vitro research use only. Not for use in diagnostics or in humans.

Distributor



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer



HOKUDO CO.,LTD

4-28,West10,Hachiken9jyo,Nishi-ku,
Sapporo,063-0849,JAPAN
TEL : +81-11-641-7507
FAX : +81-11-644-9209
URL : <http://www.hokudo.co.jp>
E-mail : info@hokudo.co.jp