



PRODUCT DATA SHEET

Product: Anti-Pyruvate dehydrogenase phosphatase 1

Cat.No: AB 41 (100 μ l , 1mg/ml)

Host:

Mouse IgG1

Species Reactivity

Rat , Mouse ,others not tested.

Immunogen

Recombinant Rat pyruvate dehydrogenase phosphatase 1 is full length amino acid and 6xHis.

Format:

100ul of purified mouse monoclonal antibody by proteinA .

Buffer : PBS(-)

Storage and Stability:

Stable for 2years when stored at less than -20°C .

Aliquot to avoid freeze/thaw cycles.

Application and Suggested Dilutions:

A working dilution of 1 : 3000 is suggested for western blots and ELISAs.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For in vitro research use only. Not for use in diagnostics or in humans.

Distributor



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer



HOKUDO CO.,LTD

4-28,West10,Hachiken9jyo,Nishi-ku,
Sapporo,063-0849,JAPAN
TEL. : +81-11-641-7507
FAX. : +81-11-644-9209
URL : <http://www.hokudo.co.jp>
E-mail : info@hokudo.co.jp