



Product Data Sheet

Anti Human HLA-A, B, C Antigen, Mouse Monoclonal antibody

- Catalog No. : HKD-AB46-EX
- Quantity : 100 μ L / vial (1mg / mL)
- Clone : EMR8-5
- Isotype : IgG1
- Buffer : Phosphate buffered saline pH7.4 , 0.05% Sodium Azide
- Immunogen : recombinant HLA-A*2402 extracellular domain
- Specificity : Heavy chains of HLA-A, HLA-B, HLA-C
- Applications : immunohistochemistry (formalin-fixed tissue , 1:50~1:100)
Western blotting
- Storage : Store at 2~8°C, For long-term storage aliquot and store at -20°C.
Avoid freeze / thaw cycles.
- References :
1. Torigoe T, Asanuma H, Shimosawa K, et al. Immunohistochemical analysis of HLA class I expression in tumor tissues revealed unusually high frequency of down-regulation in breast cancer tissues. J. Translational Med. Submitted
2. Komori H, Nakatsura T, Senju S, Yoshitake Y, Motomura Y, Beppu T, Matsui M, Torigoe T, Sato N, Baba H, Nishimura Y. Identification of HLA-A2 or -A24-restricted CTL epitopes possibly useful for glypican-3-specific immunotherapy for hepatocellular carcinoma. Clinical Cancer Res. in press

For research and in vitro use only. Not for diagnostic or therapeutic work.

Distributor



COSMO BIO CO., LTD.
Inspiration for Life Science
TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer



HOKUDO CO.,LTD
4-28,West10,Hachiken9jyo,Nishi-ku,
Sapporo,063-0849,JAPAN
TEL : +81-11-641-7507
FAX : +81-11-644-9209
URL : <http://www.hokudo.co.jp>
E-mail : info@hokudo.co.jp