




Code No.KAL-KE094-EX

For research use only

## Anti Human Urate Transporter 1 (URAT1) Polyclonal Antibody

Package Size	25 $\mu$ g (250 $\mu$ L / vial)
Format	Rabbit polyclonal antibody 0.1mg/mL
Buffer	PBS [containing 50% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Store below -20°C Once thawed, store at 4°C. Repeated freeze-thaw cycles should be avoided.
Purification method	This antibody was purified from rabbit serum by affinity chromatography.
Working dilution	For Immunohistochemistry ; 1~5 $\mu$ g/mL For Western blotting ; 0.1~0.5 $\mu$ g/mL
HGNC* Name	SLC22A12 (solute carrier superfamily 22A12 )

\*HGNC:Human Gene Nomenclature Committee



### Immunohistochemistry

Sample: Human proximal tubule

Preparation of antibodies and instruction:

Kanai Y, Anzai N

Department of Pharmacology and Toxicology,  
Kyorin University School of Medicine



Code No.KAL-KE094-EX

## Anti Human Urate Transporter 1 (URAT1) Polyclonal Antibody

### 【Reference】

1. Anzai N., Miyazaki H., Noshiro R., Khamdang S., Chairoungdua A., Shin HJ., Enomoto A., Sakmoto S., Hirata T., Tomita K., Kanai Y., Endou H. : The Multivalent PDZ Domain-Containing Protein PDZK1 Regulates Transport Activity of Renal Urate-Anion Exchanger URAT1 via Its C-terminal *J Biol Chem.*2004 Aug 10 in press
2. Ichida K., Hosoyamada M., Hisatome I., Enomoto A., Hikita M., Endou H., Hosoya T. : Clinical and molecular analysis of patients with renal hypouricemia in Japan-influence of URAT1 gene on urinary urate excretion. *J Am Soc Nephrol.*Jan;15(1).164-73, 2004
3. Koepsell H., Endou H. : The SLC22 drug transporter family. *Pflugers Arch.*Feb;447(5).666-76, 2004. Epub 2003 Jul 19.
4. Anzai N., Miyazaki H., Endou H. : Molecular biology of urate transporter. *Seikagaku* Feb;76(2).101-10, 2004 Review. Japanese. No abstract available
5. Tanaka M., Itoh K., Matsushita K., Wakita N., Adachi M., Nonoguchi H., Kitamura K., Hosoyamada M., Endou H., Tomita K. : Two male sibilings with hereditary renal hypouricemia and exercise-induced ARF. *Am J Kidney Dis.*Dec;42(6).1287-92, 2003
6. Enomoto A., Niwa T., Kanai Y., Endou H. : Urate transporter and renal hypouricemia. *Rinsyo Byori.* Sep;51(9).892-7.Japanese, 2003
7. Enomoto A., Kimura H., Chairoungdua A., shigeta Y., Jutabha P., Cha SH., Hosoyamada M., Takeda M., Skine T., Igarashi T., Mtsuo H., Kikuchi Y., oda T., Ichida K., Hosoya T., Shimokata K., Niwa T., Kanai Y., Endou H. : Molecular identification of a renal urate anion exchanger that regulates blood urate levels. *Nature.*May 23;417(6887).447-52, 2002. Epub 2002 Apr 14.

#### Distributor



COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME  
KOTO-KU, TOKYO 135-0016, JAPAN  
TEL : +81-3-5632-9617  
FAX : +81-3-5632-9618  
URL : <http://www.cosmobio.co.jp/>  
e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

#### Manufacturer



TransGenic Inc.

7-1-6 Minatojimaminami-machi,  
Chuo-ku, Kobe,650-0047 JAPAN  
TEL : +81-78-306-0590  
FAX : +81-78-306-0589  
URL : <http://www.transgenic.co.jp/>  
e-mail : [techstaff@transgenic.co.jp](mailto:techstaff@transgenic.co.jp)