



Anti PEDF Monoclonal Antibody

(Clone No. 14K)

【Reference】

1. Stratikos, E. et al. :Recombinant human pigment epithelium-derived factor (PEDF): Characterization of PEDF overexpressed and secreted by eukaryotic cells. *Protein Science* 5: 2575-2582, 1996.
2. D. W. Dawson. et al. :Pigment epithelium-derived factor: a potent inhibitor of angiogenesis. *Science* 285: 245-248, 1999.
3. Joachim Spranger, Martin Osterhoff, et al.: Loss of the antiangiogenic pigment epithelium-derived factor in patients with angiogenic eye disease. *Diabetes*.50: 2641-2645, 2001.

* These references are the background of PEDF , and are not this antibody examples .

Distributor



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer

 Trans Genic Inc.

7-1-6 Minatojimmaminami-machi,
Chuo-ku, Kobe,650-0047 JAPAN
TEL : +81-78-306-0590
FAX : +81-78-306-0589
URL : <http://www.transgenic.co.jp/>
e-mail : techstaff@transgenic.co.jp