



Code No.KAL-KR059-EX


For Research Use Only

Anti PDX1 Polyclonal Antibody, Rabbit

Pancreas duodenum homeobox 1 (PDX1) is a transcription factor encoded by homeodomain gene. This protein promotes the transcription of the pancreatic genes such as insulin and somatostatin, so that it is indispensable for keeping normal function of pancreatic beta cell. As it is involved in secretion of insulin, PDX1 becomes therapeutic target against pancreatic diseases such as diabetes and cancer.

This antibody specifically reacts with PDX1 in mouse pancreas and foregut (primordial pancreas). This antibody is useful for immunohistochemistry and could become good research tool for studying pancreatic function, embryology, and diabetes.

Package Size	25 μ g (100 μ L/vial)
Format	Rabbit polyclonal antibody 0.25mg/mL
Buffer	Block Ace as a stabilizer, containing 0.1% Proclin as a bacteriostat
Storage	Below -20°C
Purification method	This antibody was purified from rabbit serum by affinity chromatography.
Working dilution	For Immunohistochemistry ; 1 \sim 5 μ g/mL





Anti PDX1 Polyclonal Antibody, Rabbit

【Reference】

1. Eur J Endocrinol 2002 Feb;146(2):129-41
2. Genes Dev 2001 15:111-27

Distributor



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer

 Trans Genic Inc.

7-1-6 Minatojimmaminami-machi,
Chuo-ku, Kobe,650-0047 JAPAN
TEL : +81-78-306-0590
FAX : +81-78-306-0589
URL : <http://www.transgenic.co.jp/>
e-mail : techstaff@transgenic.co.jp