




Code No.KAL-KR064-EX

For research use only

## Anti Rat Podocalyxin Polyclonal Antibody, Rabbit

Package Size	25 $\mu$ g (100 $\mu$ L/vial)
Format	Rabbit polyclonal antibody 0.25mg/mL
Buffer	PBS [containing 50% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Below -20 $^{\circ}$ C
Purification method	This antibody was purified from rabbit serum immunized with synthesized partial peptide of Rat Podocalyxin by peptide affinity chromatography.
Working dilution	Western Blotting : 2~5 $\mu$ g/mL Immunofluorescence : 10 $\mu$ g/mL



### Western blotting

Sample : 200 $\mu$ g/lane

- ① Rat liver lysate (Negative control)
- ② Rat kidney lysate
- ③ Rat heart lysate

Preparation of antibodies and instruction :

Dr. Takeda T. and Dr. Gejyo F. at Department of Applied Molecular Medicine, Niigata University Graduate School of Medical and Dental Sciences



## Anti Rat Podocalyxin Polyclonal Antibody, Rabbit

### 【Reference】

1. Takeda T. et al. : J Clin Invest 2001 Jul;108(2):289-301
2. Kerjaschki D. et al. : Am J Pathol 1985 Mar;118(3):343-9
3. Kerjaschki D. et al. : J Cell Biol 1984;98:1591-6

### Distributor



COSMO BIO Co., LTD.  
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME  
KOTO-KU, TOKYO 135-0016, JAPAN  
TEL : +81-3-5632-9617  
FAX : +81-3-5632-9618  
URL : <http://www.cosmobio.co.jp/>  
e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

### Manufacturer

 Trans Genic Inc.

7-1-6 Minatojiminami-machi,  
Chuo-ku, Kobe, 650-0047 JAPAN  
TEL : +81-78-306-0590  
FAX : +81-78-306-0589  
URL : <http://www.transgenic.co.jp/>  
e-mail : [techstaff@transgenic.co.jp](mailto:techstaff@transgenic.co.jp)