





Code No.KAL-KR074-EX

For research use only

Anti Keratin 12 Polyclonal Antibody

Package Size	25 μ g (100 μ L/vial)
Format	Rabbit polyclonal antibody 0.25mg/mL
Buffer	PBS [containing 50% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Below -20°C
Purification method	This antibody was purified from rabbit serum immunized with synthesized partial peptide of rabbit Keratin 12 by peptide affinity chromatography.
Working dilution	Immunohistochemistry : 2 μ g/mL



Immunohistochemistry

Sample : Human Cornea

Preparation of antibodies and instruction :

Dr. Tano Y., Dr. Nishida K. and Dr. Watanabe K.

at Department of Ophthalmology, Osaka University Medical School



Anti Keratin 12 Polyclonal Antibody

【Reference】

1. Liu CY. et al. : Curr Eye Res. 1993 Nov;12(11):963-74
2. Nishida K. et al. : Invest Ophthalmol Vis Sci. 1996 Aug;37(9):1800-9
3. Moyer PD. et al. : Differentiation. 1996 Mar;60(1):31-8.
4. Kinoshita S. et al. : Prog Retin Eye Res. 2001 Sep;20(5):639-73

Distributor



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer

 Trans Genic Inc.

7-1-6 Minatojimaminami-machi,
Chuo-ku, Kobe,650-0047 JAPAN
TEL : +81-78-306-0590
FAX : +81-78-306-0589
URL : <http://www.transgenic.co.jp/>
e-mail : techstaff@transgenic.co.jp