




Code No.KAL-KR083-EX

For research use only

Anti Matrix Gla Protein Polyclonal Antibody

Gla : gamma-carboxyglutamic acid

Package Size	50µg (200µL/vial)
Format	Rabbit polyclonal antibody 0.25mg/mL
Buffer	PBS [containing 50% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Below -20°C
Purification method	This antibody was purified from rabbit serum immunized with synthesized partial peptide of mouse Matrix Gla Protein by peptide affinity chromatography.
Working dilution	Immunohistochemistry : 15µg/mL



Immunohistochemistry

Sample : Mouse Limb (Embryo 14 days)
(Calcification site and its around mature hypertrophic chondrocyte are positively stained.)

Preparation of antibodies and instruction

Dr. Mano H. and Dr. Nakatani S. at Department of Clinical Dietetics and Human Nutrition, Faculty of Pharmaceutical Sciences, Josai University



Anti Matrix Gla Protein Polyclonal Antibody

Gla : gamma-carboxyglutamic acid

【Reference】

1. Canfield AE. et al. : FEBS Lett. 2000 Dec 29;487(2):267-71
2. Wallin R. et al. : Med Res Rev. 2001 Jul;21(4):274-301

Distributor



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer

 Trans Genic Inc.

7-1-6 Minatojimaminami-machi,
Chuo-ku, Kobe, 650-0047 JAPAN
TEL : +81-78-306-0590
FAX : +81-78-306-0589
URL : <http://www.transgenic.co.jp/>
e-mail : techstaff@transgenic.co.jp