




Code No.KAL-KR099-EX

For research use only

Anti Human 17 β HSD 8 Polyclonal Antibody

17 β HSD 8 : 17 β hydroxysteroid dehydrogenase type 8

| | |
|---------------------|--|
| Package Size | 25 μ g (100 μ L / vial) |
| Format | Rabbit polyclonal antibody 0.25mg/mL |
| Buffer | PBS [containing 50% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat] |
| Storage | Store below -20 $^{\circ}$ C Once thawed, store at 4 $^{\circ}$ C. Repeated freeze-thaw cycles should be avoided. |
| Purification method | This antibody was purified from rabbit serum immunized with synthesized partial peptide of human 17 β HSD8 by peptide affinity chromatography. |
| Working dilution | For Immunohistochemistry ; 5 μ g/mL |



Immunohistochemistry

Sample : Human endometrium

Preparation of antibodies and instruction :

Dr. Okamura H. at Department of Reproductive Medicine and Surgery, Faculty of Medical and Pharmaceutical Sciences Kumamoto University



Code No.KAL-KR099-EX

Anti Human 17 β HSD 8 Polyclonal Antibody

17 β HSD 8 : 17 β hydroxysteroid dehydrogenase type 8

【Reference】

1. Ramirez S. et al. : Mol Cell Endocrinol. 1998 Aug 25;143(1-2):9-22
2. Woo D. et al. : Mol Cell Endocrinol. 2001 May 15;176(1-2):155-62

Distributor



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer



TransGenic Inc.

7-1-6 Minatojiminami-machi,
Chuo-ku, Kobe,650-0047 JAPAN
TEL : +81-78-306-0590
FAX : +81-78-306-0589
URL : <http://www.transgenic.co.jp/>
e-mail : techstaff@transgenic.co.jp