




Code No.KAL-KR100-EX

For research use only

## Anti Mouse Semaphorin 3B Polyclonal Antibody

Package Size	100 $\mu$ g (400 $\mu$ L / vial)
Format	Rabbit Polyclonal antibody 0.25mg/mL
Buffer	PBS [containing 50% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Store below -20 $^{\circ}$ C Once thawed, store at 4 $^{\circ}$ C. Repeated freeze-thaw cycles should be avoided.
Purification method	This antibody was purified from rabbit serum immunized with synthesized partial peptide of Mouse Semaphorin 3B by peptide affinity chromatography.
Use	Immunohistochemistry, Immunoprecipitation, Western blotting
Working dilution	For Immunohistochemistry ; 15 $\mu$ g/mL For Western blotting ; 8 $\mu$ g/mL For Immunoprecipitation ; 2 $\mu$ g



### Immunohistochemistry

Sample : Human placenta

Preparation of antibodies and instruction :

Dr. Suzuki, T. at Second Department of Pathology, School of  
Medicine, Fukushima Medical University



Code No.KAL-KR100-EX

## Anti Mouse Semaphorin 3B Polyclonal Antibody

### 【Reference】

1. Castro-Rivera E. et al. : Proc Natl Acad Sci USA. 2004 Aug 3;101(31):11432-7
2. Tse C. et al. : Cancer Res. 2002 Jan 15;62(2):542-6
3. Gong S. et al. : Nature. 2003 Oct 30;425(6961):917-25

#### Distributor



COSMO BIO Co., LTD.  
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME  
KOTO-KU, TOKYO 135-0016, JAPAN  
TEL : +81-3-5632-9617  
FAX : +81-3-5632-9618  
URL : <http://www.cosmobio.co.jp/>  
e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

#### Manufacturer



Trans Genic Inc.

7-1-6 Minatojimaminami-machi,  
Chuo-ku, Kobe,650-0047 JAPAN  
TEL : +81-78-306-0590  
FAX : +81-78-306-0589  
URL : <http://www.transgenic.co.jp/>  
e-mail : [techstaff@transgenic.co.jp](mailto:techstaff@transgenic.co.jp)