




Code No.KAL-KR111-EX

For research use only

Anti human Tensin2 Polyclonal Antibody

Package Size	50 μ g (200 μ L / vial)
Format	Rabbit polyclonal antibody 0.25mg/mL
Buffer	PBS [containing 50% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Store below -20 $^{\circ}$ C Once thawed, store at 4 $^{\circ}$ C. Repeated freeze-thaw cycles should be avoided.
Purification method	This antibody was purified from rabbit serum immunized with synthesized partial peptide of human Tensin2 by peptide affinity chromatography.
Working dilution	For Immunohistochemistry ; 5 μ g/mL For Western blotting ; 2 μ g



Immunohistochemistry

Sample : Mouse Kidney

Preparation of antibodies and instruction :

Dr.Agui T. and Sasaki N.

The Graduate School of Veterinary Medicine, Hokkaido University Department of Disease Control



Code No.KAL-KR111-EX

Anti human Tensin2 Polyclonal Antibody

【Reference】

1. Lo SH. : Int J Biochem Cell Biol. 2004 Jan;36(1):31-4.
2. Lo SH. et al. : J Cell Biol. 1997 Mar 24;136(6):1349-61.
3. Calderwood DA et al. : Proc Natl Acad Sci U S A. 2003 Mar 4;100(5):2272-7.
4. Chen H. : Proc Natl Acad Sci U S A. 2002 Jan 22;99(2):733-8.
5. Lo SH. : Cancer Res. 2002 Aug 1;62(15):4217-21.

Distributor



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer



TransGenic Inc.

7-1-6 Minatojimaminami-machi,
Chuo-ku, Kobe,650-0047 JAPAN
TEL : +81-78-306-0590
FAX : +81-78-306-0589
URL : <http://www.transgenic.co.jp/>
e-mail : techstaff@transgenic.co.jp