



Product Information

CAS Registry No. 327-97-9

Product Name: Chlorogenic acid

Catalog No. C 2944

Chemical Name: [1S-(1 α ,3 β ,4 α ,5 α)]-3-[[3-(3,4-Dihydroxyphenyl)-1-oxo-2-propenyloxy]-1,4,5-trihydroxycyclohexane-carboxylic acid

Synonym(s): 3-Caffeoylquinic acid

Mol. formula: C₁₆H₁₈O₉

Formula wt: 354.31

M.p. 207°C-209°C (dec.)

B.p.

Purity: ≥98%


Solubility: Soluble in ethanol or acetone.

Stability (in original shipping container):

Storage temp: Ambient

Product Description:

A chemopreventive analog of caffeic acid.



References:

Tanaka, T., Kojima, T., Kawamori, T., et al. Carcinogenesis. 4:1321-1325 (1993).

PACKING

100 mg

250 mg

1 g

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution:

This product is intended for laboratory and research use only. It is not for human or drug use.