



COSMO BIO CO., LTD.  
Inspiration for Life Science

## MOUSE ANTI Ca19.9 MONOCLONAL ANTIBODY

|                         |                                                                                                                                                        |
|-------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No.</b>      | LNM-KR-014                                                                                                                                             |
| <b>Clone</b>            | 3C11                                                                                                                                                   |
| <b>Concentration</b>    | 1 mg/ml                                                                                                                                                |
| <b>Volume</b>           | 100 UL                                                                                                                                                 |
| <b>Subclass</b>         | IgA                                                                                                                                                    |
| <b>Immunogen</b>        | Human CA19.9                                                                                                                                           |
| <b>Specificity</b>      | Purified CA19.9 from tumor cells                                                                                                                       |
| <b>Applications</b>     | Sandwich ELISA( This one clone is useful for<br>both labeled antibody and capture antibody.)                                                           |
| <b>Format</b>           | Purified immunoglobulin from BALB/C Mouse Ascites<br>(Salting -out fractionation by ammonium sulfate and<br>ion exchange chromatographic purification) |
| <b>Formulation</b>      | Liquid in 15mM Phosphate buffer containing 0.15M NaCl,<br>0.09% sodium azide, pH 7.2                                                                   |
| <b>Storage/Handling</b> | 2 to 8°C, Maintain at -70°C for up to 12 months.<br><br>Avoid repeated freeze/thaw cycle.                                                              |

*For research use only; not for use as a diagnostic*



COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN  
TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp

[www.cosmobio.co.jp](http://www.cosmobio.co.jp)

## 抗ヒトCa19.9モノクローナル抗体

|       |                                                              |
|-------|--------------------------------------------------------------|
| 品 番   | K R - 0 1 4                                                  |
| クローン  | 3 C 1 1                                                      |
| 型     | I g A                                                        |
| 由 来   | B A L B / C マウス腹水                                            |
| 精 製 法 | 硫酸塩析, r - プロテインA                                             |
| 溶 液   | P B S (0.09%アジ化ナトリウム含む)                                      |
| 容 量   | 1 0 0 $\mu$ g                                                |
| 保 存 法 | 4 ~ 8℃<br>長期保存の場合は、- 7 0℃以下にて凍結保存して下さい。<br>凍結融解は出来るだけ避けて下さい。 |
| 備 考   | 抗ヒトCa19.9抗体(3C11)はサンドイッチELISA法に適しています。                       |



コスモ・バイオ株式会社