

抗ヒトC a 19.9 モノクローナル抗体

品 番	K R - 0 1 4
クローン	3 C 1 1
型	I g A
由 来	培養上清
精 製 法	クロマトグラフィー
溶 液	P B S (0.09%アジ化ナトリウム含む)
数 量	1 0 0 μ g
保 存 法	4 ~ 8℃ 凍結融解は出来るだけ避けて下さい。
備 考	1 クローンでサンドイッチ ELISA が可能です。



人と科学のステキな未来へ

コスモ・バイオ株式会社

— 商品に関するお問い合わせ —

TEL: 03-5632-9610 FAX: 03-5632-9619

本社所在地 〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル

**MOUSE ANTI-CA19.9
MONOCLONAL ANTIBODY**

Catalog No.	LNМ-KR-014
CLONE:	3C11
QUANTITY:	0.1mg
SPECIFICITY:	Human CA19.9
IMMUNOGEN:	Purified CA19.9 from tumor cells
Type/Subtype:	IgA
APPLICATIONS:	Sandwich ELISA(This one clone is useful for both labeled antibody and capture antibody.)
FORMAT:	Purified immunoglobulin from culture supernatant (Ion exchange chromatographic purification)
PRESENTATION:	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
STORAGE/HANDLING:	2 to 8°C, Avoid repeated freeze/thaw cycle.

For research use only; not for use as a diagnostic.



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN
TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618 / e-mail: export@cosmobio.co.jp / URL: www.cosmobio.com