



MOUSE ANTI CALPROTECTIN MONOCLONAL ANTIBODY

Catalog No.	LNM-KR-015
Clone	No.134
Size	0.1mg
Subclass	IgG1,k
Immunogen	Human Calprotectin
Specificity	Human Calprotectin
Applications	Sandwich ELISA (This antibody can be used as a capture antibody in a human Calprotectin ELISA in combination with biotinlated, monoclonal antibody No.19.)
Format	Purified immunoglobulin from BALB/C Mouse Ascites Anmmonium sulfate fractionation HPLC
Formulation	PBS: Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
Storage/Handling	4 to 8°C, Maintain at -70°C for up to 12 months. Avoid repeated freeze/thaw cycle.
Note	IgG concentration of each lot is described in the URL below. https://www.cosmobio.com/products/lnm_kr015_01.html

For research use only; not for use as a diagnostic



抗ヒトカルプロテクチンモノクローナル抗体

品 番	K R - 0 1 5
クローン	N o . 1 3 4
型	I g G 1 , κ
由 来	B A L B / C マウス腹水
精 製 法	硫酸塩析, r - プロテインA
溶 液	P B S (0.09%アジ化ナトリウム含む)
容 量	1 0 0 μ g
保 存 法	4 ~ 8 °C 長期保存の場合は、- 7 0 °C以下にて凍結保存して下さい。 凍結融解は出来るだけ避けて下さい。
備 考	抗ヒトカルプロテクチン抗体 (No.134) はサンドイッチ ELISA 法の固相に適しています。なお、二次抗体としては抗ヒトカルプロテクチン抗体 (No.19) が適しています。



コスモ・バイオ株式会社