

## 抗ヒトミオグロビンモノクローナル抗体

品 番	KR-028
クローン	No. 44
型	IgG1, $\kappa$
由 来	BALB/Cマウス腹水
精製法	硫酸塩析, r-プロテインA
溶 液	PBS (0.09%アジ化ナトリウム含む)
数 量	100 $\mu$ g
保 存 法	4~8℃ 凍結融解は出来るだけ避けて下さい。
備 考	抗ヒトミオグロビン抗体 (No.44) はサンドイッチ ELISA 法の二次抗体に適しています。なお、固相としては抗ヒトミオグロビン抗体 (No.29) が適しています。



人と科学のステキな未来へ

コスモ・バイオ株式会社

— 商品に関するお問い合わせ —

TEL: 03-5632-9610 FAX: 03-5632-9619

本社所在地 〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル

**MOUSE ANTI MYOGLOBIN  
MONOCLONAL ANTIBODY**

**Catalog No.** LNM-KR-028  
**CLONE:** No.44  
**QUANTITY:** 0.1mg  
**SPECIFICITY:** Human Myoglobin  
**IMMUNOGEN:** Human Myoglobin  
**Type/Subtype:** IgG1  $\kappa$   
**APPLICATIONS:** Sandwich ELISA (This antibody can be used as a detection antibody in a human Myoglobin ELISA combination with the monoclonal capture antibody No.29.)  
**FORMAT:** Purified immunoglobulin from BALB/C Mouse Ascites (Affinity chromatography)  
**PRESENTATION:** Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2  
**STORAGE/HANDLING:** 2 to 8°C, Avoid repeated freeze/thaw cycle.

***For research use only; not for use as a diagnostic.***



**COSMO BIO CO., LTD.**  
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN  
TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618 / e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp) / URL: [www.cosmobio.com](http://www.cosmobio.com)